





4K Ultra HD Endoscope

Highest brilliance for excellent results

- Uniform brightness and visibility in the entire field of view
- Higher precision for more safety
- Sapphire glass distal and proximal
- Ø 4 mm, 175 mm WL



3840 x 2160

Superb sharpness of the finest texture due to 4x higher resolution



1920 x 1080

SPIGGLE & THEIS
Medizintechnik GmbH

Burghof 14
51491 Overath / Germany
Tel.: +49 (0) 2206 / 90 81 - 0
Fax: +49 (0) 2206 / 90 81 -13
info@spiggle-theis.de
www.spiggle-theis.de



FL_4K_EN_01



PO Box 3403, Erina NSW 2250, Australia
P +61 2 9810 6947 **F** +61 2 8072 1556
admin@coremed.com.au **www.coremed.com.au**